

## **SOP - Extraction of non-volatile, small, polar metabolites**

The entire 400-ml cultures were centrifuged ( $17,700 \times g$  for 15 min at  $4^{\circ}\text{C}$ ). Cell pellets were immediately washed with 100 ml NaCl (0.9 % [wt/vol], anoxic for anaerobic cultures) and centrifuged ( $17,700 \times g$  for 15 min at  $4^{\circ}\text{C}$ ). The resultant pellets were immediately resuspended in 2 ml MeOH containing 0.02–0.05 mM ribitol, transferred in portions á 1 ml into Precellys<sup>®</sup> tubes (see above), shock frozen in liquid  $\text{N}_2$  and stored at  $-80^{\circ}\text{C}$  until further analysis.

Samples were thawed on ice before homogenization at  $-10^{\circ}\text{C}$  with 3 cycles at a speed of 6,800 rpm for 30 s with breaks of 30 s, using Precellys24 homogenizer (Bertin Instruments, Montigny-le-Bretonneux, France). Afterwards, 400  $\mu\text{l}$  water were added to the sample, followed by a mixing step of 2 min at 2,000 rpm (MixMate; Eppendorf, Hamburg, Germany). Then, 250  $\mu\text{l}$  chloroform were added and the sample mixed again for 2 min at 2,000 rpm prior to centrifugation (10,000 rpm for 10 min at  $4^{\circ}\text{C}$ ). Subsequently, the polar phase was then transferred to 2-ml reaction tubes and dried at  $14^{\circ}\text{C}$  under vacuum. Samples were stored at  $-20^{\circ}\text{C}$  until derivatization for GC-EI-MS analysis. For every growth condition 4 biological samples with 2 technical replicates each were analyzed.